

icare
HOME



IOP SELF-MEASURING
ANYWHERE, ANYTIME



NO ANESTHESIA
NO AIR PUFF



THE 24HR TONOMETER

The easy-to-use Icare® tonometers revolutionize IOP measurements making it quick, effortless and effective. The Icare® HOME tonometer is designed for patients who need additional IOP monitoring according to their doctor's recommendation.

The device is based on a rebound measuring principle that eye care professionals have come to trust. A small probe touches the cornea and the impact duration, deceleration and other motion parameters are evaluated.

Icare® HOME does not require topical anesthetic or specialized skills for use. The HOME tonometer records IOP measurements and additional data which is not displayed to the patient. The healthcare professional downloads the measurement data via Icare LINK software.

KEY FEATURES



Icare® AMS

Automatic measuring sequence:
series and single mode
with one button



Icare® EasyPos

Intelligent positioning
assistant for the
correct alignment
of the tonometer




Icare® EyeSmart

Automatic OD/OS
recognition technology

IOP SELF-TONOMOMETRY WITH ICARE HOME

- + REAL WORLD DATA
- + IDENTIFY IOP VARIATION AND SPIKES
- + VALUABLE ADJUNCT TO CLINIC

mmHg



CLINIC MEASUREMENTS (12 months)

- IOP monitoring in the clinic
- 24 hour IOP self-measuring with Icare HOME

mmHg



Icare HOME (24 HOURS)



ICARE USA - VISIONARY IN VISION

Icare USA is a subsidiary of Icare Finland, the original developer and global leader in handheld tonometry. Our patented technology (more than 20 patents and patent applications) combined with ISO 13485 certified quality system has made us a respected player in our field of expertise.

The advanced Icare® product line offers reliable, precise and reproducible accuracy when measuring intraocular pressure in both experienced and inexperienced hands. Icare tonometers available:

Icare® HOME, Icare® ic100, Icare® TA01i

Made in Finland.

Taking control of your eye health is as easy as 1,2,3!

"...Renting this device has made a huge difference in managing my glaucoma now that I know my IOP fluctuates so much..."

1. Obtain a prescription from your Doctor
2. Visit MyEYES.net
3. Fill out an online order form & Rent (or purchase) a device

We will ship it directly to your door!

Simple, Safe & Secure

Here at MyEYES, we are dedicated to assisting you with your eye health. We care and understand the difficulties of glaucoma. That's why using this device can help you literally take back control of your eye health. IOP can be a vital link in your treatment protocols and by using this device, you can take CONTROL!

ICARE HOME APP



Find us at:

MyEYES.net

Rx Only